

BioMasher II World First Disposable Homogenizer

Feature

- Dimple Processing

The stir stick and inside of micro tube are dimple processing, so it is best homogenizing efficiency. (The homogenizing efficiency is about twice compared with past one)

- Repercussion Prevention

The repercussion prevention sleeve is on the stir stick, so the sample is prevented from rebounding besides the micro tube.

Application

- Homogenizing of animal tissue and organs
- Homogenizing of leaf, stalk, root and seed of plant
- Researches such as insect, fish, shellfish, bacillus, yeasts, DNA, RNA, protein and enzyme etc.

Material



1.5ml Micro Tube: PP (Polypropylene)

Stir Stick: POM (Polyacetal)

Ordering Number

| Item No. | Description | Q'ty |
|----------|----------------------------------|------------|
| 890863 | Biomasher II , Sterilization* | 50 sets/pk |
| 890853 | Biomasher II , Non-Sterilization | 50 sets/pk |

* EOG Sterilized




Power Masher

For BioMasher II

Feature

- Speed Adjust: 3 steps
- Water Proof Type
- Battery: AA type x 1



Item no.: 320105